

Test With Changes in Collection - Quest 11/19/25 Go-live

Below is a crosswalk of the testing scheduled to go-live on 11/19/2025 that have changes to their collection containers. Links to the Quest lab test directory are provided below. Details on these tests will be updated in our lab test directories on 11/19.

[Laboratory Catalogs | Corewell Health](#)

Current EPIC Code	FUTURE Quest EPIC Order	CPT	Collection Detail	Preferred Collection	Container Changes
CREATINE, SERUM [LAB1231186]	CREATINE, SERUM [LAB1231958]	82540	6 ML RED TOP F [1230102021]	Only frozen serum from Red-top	Plasma no longer acceptable
CYP2C19 GENOTYPING [LAB1231231]	CYP2C19 GENOTYPE [LAB1231987]	81225	4 ML LAVENDER TOP R [1230300351]	Lav/Pink EDTA whole blood	No longer ACD
DEXAMETHASONE, SERUM [LAB1230895]	DEXAMETHASONE [LAB1231936]	80299	6 ML RED TOP R [1230300000]	Red top only	No longer SST
INTERLEUKIN 2 RECEPTOR (CD25), SOLUBLE [LAB1230002]	INTERLEUKIN-2 RECEPTOR ALPHA CHAIN (IL-2RA/CD25), SOLUBLE [LAB1231942]	83520	6 ML RED TOP F [1230102021]	Only frozen serum from Red-top	No longer SST
NUCLEAR MAGNETIC RESONANCE LIPOPROTEIN PROFILE, SERUM [LAB1230522]	LIPOPROTEIN FRACTIONATION ION MOBILITY [LAB1231939]	83704	5 ML GOLD TOP SST R [1230300004]	2x Red top tube	No longer required to use NMR kit; No longer SST
PM 1 ANTIBODY [LAB1231217]	PM/SCL-100 ANTIBODY, IGG BY IMMUNOBLOT [LAB1231989]	86235	5 ML GOLD TOP SST R [1230300004]	SST tube	No longer Red-top tube
PORPHOBILINOGEN, URINE 24 HR [LAB1231235]	PORPHOBILINOGEN, QUANTITATIVE, 24-HOUR URINE [LAB1231991]	84110	24 HR URINE CONTAINER R [1230300322]	24HR urine container refrigerated	No preservative needed
THIOPURINE METHYLTRANSFERASE (TPMT) ENZYME LEVEL [LAB10144]	TPMT ACTIVITY [LAB1232000]	84433	4 ML LAVENDER TOP R [1230300351] x2	2x Lavender EDTA whole blood	No longer Dark Green Top
VITAMIN K1 [LAB1231467] CHE only to Mayo	VITAMIN K [LAB3334]	84597	3 ML LAVENDER TOP F [1230102002]	Lavender EDTA Plasma - Protect from Light / Frozen	No longer required to use Mayo Lavender collection tube